










Technical Manual



MIAMI LUMBAR[®] with Thoracic Extensions

		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	18
IT	Istruzioni per l'uso	22
NO	Bruksanvisning	26
DA	Brugsanvisning	30
SV	Bruksanvisning	34
FI	Käyttöohjeet	38
NL	Gebruiksaanwijzing	42
PT	Instruções de Utilização	46
JA	取扱説明書	50
ZH	中文说明书	53

LIST OF PARTS

- a. Back Panel with Compression System
- b. Compression System Handles
- c. Overlap Closures
- d. Yellow Attachment Strap
- e. Blue Attachment Strap
- f. Posterior Panel with Lateral Support
- g. Posterior Panel without Lateral Support
- h. Axillary Strap Guide
- i. Axillary Strap
- j. Axillary Strap Padding
- k. Posterior Thoracic Extension (PTE)
- l. PTE Strap Holder
- m. Rigid Anterior Panel with Anterior Thoracic Extension (ATE)
- n. Anterior Panel Slot
- o. Quick-Release Buckle
- p. Pectoral Pads
- q. Hinged Adjustment Mechanism
- r. Black Allen Wrench
- s. Black Posts
- t. Silver Allen Wrench
- u. PTE Screws
- v. PTE Plugs

INDICATIONS FOR USE

Spinal stenosis, Herniated disc, Degenerative spinal pathologies, Post-surgical stabilization, Stable, non-displaced spinal fractures, Spondylolithesis, Spondylolysis

CAUTION

- The Miami Lumbar with anterior and/or posterior thoracic extension(s) is a supportive brace only and is not intended or guaranteed to prevent spinal injury. It should be applied by a licensed healthcare practitioner and as directed by a physician or other qualified medical authority. Ossur accepts no responsibility for any injury which might occur while using the orthosis.
- Do not over tighten the compression system to the point where it causes discomfort or difficulty breathing. For sitting, the patient may find it desirable to slightly loosen the compression system.
- If brace causes discomfort or irritation, contact your healthcare provider for further instructions.

SIZE SELECTION

The thoracic attachments are size adjustable and compatible with all sizes of the Miami Lumbar. Please see the Miami Lumbar Instructions for Use to determine the correct size belt to be used in your application.

BRACE ASSEMBLY

1. Size the Miami Lumbar belt according to the Miami Lumbar Instructions for Use.
2. In advance of applying the PTE, remove the top portion of the padded cover from the posterior panel to expose the two holes labeled “PTE”. Starting from the outside of the panel, insert the two black posts and secure with the PTE screws (applied from the inside of the panel).
(NOTE: The screws may have shipped already inserted in the post. Remove screw first if necessary). Tighten with the provided silver allen wrench (**Fig. 1**).
3. To attach the posterior panel to the belt, thread the yellow attachment strap through the “H” shaped slot on the posterior panel (**Fig. 2**). For different panel-belt attachment types and methods, please refer to QRC 0038, supplied as an insert along with this manual.
4. The TLSO product is supplied with the Pectoral/Axillary Strap configuration (**Fig. 3**). Other brace configurations

are possible with this product:

- Over-The-Shoulder Strap configuration (**Fig. 4**): Remove the axillary straps from their keyholes on the shoulder portion of the PTE and discard.
- PTE only brace configuration (**Fig. 5**): This product includes the PTE and axillary straps only (does not include the shoulder straps or ATE). **note:** For application of this brace configuration, follow instructions below disregarding any steps pertaining to the ATE.

BRACE APPLICATION – SUPINE

NOTE: Össur recommends applying this product while the patient is in a supine position.

1. Logroll the patient on to their side.
2. Position the back panel centered over the spine with the bottom of the rigid posterior panel at approximately the sacrococcygeal junction (**Fig. 6**).
3. To size the PTE, position it so that the shoulder extensions rest on top of the upper trapezius. Note the number corresponding to the keyhole which lines up with the black posts in the posterior panel (**Fig. 7a**). Remove the brace from the patient and assemble the PTE by aligning the black posts with the keyholes in the PTE that correspond with the correct number and slide the PTE sideways to lock in place (**Fig. 7b**). Place the PTE plugs into the hole next to the PTE post in order to prevent unintended dislodging (**Fig. 7c**)
NOTE: for shorter patients, the PTE can be slid behind the yellow attachment strap on the lumbar belt.
4. While the patient is on their side, take the long end of the axillary strap on the upper shoulder and wrap it through the patient's axilla, adjusting the padding as necessary. Bring it across the back and tuck the end under the patient at the level of the waist. The strap should affix to the belt at this location (**Fig. 8**). Log roll the patient back into the supine position.
5. Wrap the belt around the patient's waist and temporarily secure the belt with the overlap closures (**Fig. 9**).
6. Log roll the patient to the opposite side and repeat step 4 with the other axillary strap. Do not tighten the straps at this time.
(**NOTE:** Disregard if not using axillary straps.)
Log roll the patient back into the supine position.
7. Loosen the overlap closure and attach the ATE to the belt by threading the blue attachment strap through the "H" shaped slot on the anterior panel. It may be necessary to slide the ATE up to expose the anterior panel slot (**Fig. 10**). For different panel-belt attachment types and methods, please refer to QRC 0038, supplied as an insert along with this manual.
8. Re-secure the overlap closures, adjusting the posterior and/or anterior panels as needed to achieve the correct fit and positioning (**Fig. 9**). The posterior panel should be centered on the spine with the bottom of the posterior panel at approximately the sacrococcygeal junction. The anterior panel should be centered on the abdomen with the bottom edge just above the symphysis pubis while still allowing the patient to sit comfortably.
9. Adjust the height of the ATE by pushing on the button in the anterior panel and sliding the strut to the appropriate length (**Fig. 11**). The superior edge of the pectoral pads should be located approximately one inch below the sternal notch.
10. Using the provided black allen wrench, adjust the angles on the two hinged adjustment mechanisms to accommodate sternal contours (**Fig. 12**).
NOTE: Make sure the hinged adjustment mechanisms are tightened sufficiently to lock them into position.
11. Re-adjust the overlapping closure if necessary. To tighten the compression system, slide thumbs through the holes in the two compression system handles and pull to the appropriate tightness (**Fig. 13**).
12. Connect both quick-release buckles and tighten the shoulder straps (**Fig. 14**).
13. Apply the axillary strap guides to both sides of the belt to serve as guides for the axillary straps during re-application (**Fig. 15**). Place the distal edge of the guides under the hook and loop connection of the belt arm to back panel (**Fig 15a**). The guides may be angled and/or trimmed to ensure straight positioning of the axillary straps. Figures 15b and 15c depict the right and wrong way of attaching the guides, respectively.
14. Tighten the axillary straps by pulling forward and attaching to the front of the brace (**Fig. 16**). If necessary, adjust the axillary strap padding for patient comfort and/or trim strap length.
15. Upon sitting, standing or ambulation, the brace may require fine-tuning for comfort, fit and stabilization. This may include adjusting or tightening straps, compression system, rigid panels, and/or ATE/PTE.

BRACE APPLICATION – STANDING

1. Position the back panel centered over the spine with the bottom of the rigid posterior panel at approximately the sacrococcygeal junction.
2. To size the PTE, position it so that the shoulder extensions rest on top of the upper trapezius. Note the number corresponding to the keyhole which lines up with the black posts in the posterior panel. Remove the brace from the patient and assemble the PTE by sliding the extension behind the yellow attachment strap on the lumbar belt. Align the black posts with the keyholes in the PTE that correspond with the correct number and slide the PTE sideways to lock in place.
3. Loosen the overlap closure and attach the ATE to the belt by threading the blue attachment strap through the “H” shaped slot on the anterior panel. It may be necessary to slide the ATE up to expose the anterior panel slot. For different panel-belt attachment types and methods, please refer to QRC 0038, supplied as an insert along with this manual.
4. Wrap the belt around the patient's waist and temporarily secure the belt with the overlap closure.
5. Take the long end of one of the axillary straps and wrap it through the patient's axilla, adjusting the padding as necessary. Bring it across the back and affix it to the belt just above the contralateral hip. Repeat for the other axillary strap. Do not tighten the straps at this time. (Disregard if not using axillary straps.)
6. Re-secure the overlap closure, adjusting the posterior and/or anterior panels as needed to achieve the correct fit and positioning. The posterior panel should be centered on the spine with the bottom of the posterior panel at approximately the sacrococcygeal junction. The anterior panel should be centered on the abdomen with the bottom edge just above the symphysis pubis while still allowing the patient to sit comfortably.
7. Adjust the height of the ATE by pushing on the button and sliding the strut to the appropriate length. The superior edge of the pectoral pads should be located approximately one inch below the sternal notch.
8. Using the provided black allen wrench, adjust the angles on the two hinged adjustment mechanisms to accommodate sternal contours.
NOTE: Make sure the hinged adjustment mechanisms are tightened sufficiently to lock them into position.
9. Re-adjust the overlapping closure if necessary. To tighten the compression system, slide thumbs through the holes in the two compression system handles and pull to the appropriate tightness.
10. Connect both quick-release buckles and tighten the shoulder straps.
11. Apply the axillary strap guides to both sides of the belt to serve as guides for the axillary straps during re-application. Refer to step 13. in the instructions for supine application for correct attachment of the axillary strap guides.
12. Tighten the axillary straps by pulling forward and attaching to the front of the brace. If necessary, adjust the axillary strap padding for patient comfort and/or trim strap length.
13. Upon sitting, standing or ambulation, the brace may require fine-tuning for comfort, fit and stabilization. This may include adjusting or tightening straps, compression system, rigid panels, and/or ATE/PTE.

BRACE REMOVAL

1. Detach axillary straps and loosen (if applicable). Affix them on the brace just above the hip (**Fig. 17**).
2. Detach compression system handles and loosen. Affix them just above the axillary straps.
3. Detach the quick-release buckle at the right shoulder (if applicable).
4. Detach the overlap closure and remove brace (**Fig. 18**).
5. Before reapplying the brace and to ensure a proper fit, be sure the compression system is stretched to its full width by pulling the ends of the belt apart (**Fig. 19**).

BRACE ADJUSTMENTS

- The posterior and anterior rigid panels can be modified as necessary to optimize patient fit and comfort. Simply remove the panels from their fabric sleeve and adjust the rigid panel with a heat gun (F/C) and/or trimming device.
- The Miami Lumbar system is completely modular and can be easily customized to the needs of your patient. The ATE, PTE, anterior panels, posterior panels and/or strapping configurations can be added or removed depending on the motion restriction and compression desired throughout the rehabilitation period.

BRACE MAINTENANCE

The Miami Lumbar system is designed for minimum maintenance and care. The padded covers may be removed from the rigid panels for washing.

- Hand-wash using mild soap and rinse thoroughly.
- Air dry flat.
- **DO NOT** machine wash or place in tumble dryer.

NOTE: The pectoral pads are able to be removed from the ATE (i.e. for cleaning). To remove, pull and twist the pectoral pads until the ball snaps out of the socket. To re-assemble, press the ball into the socket until it snaps back in place.

ACCESSORIES/OPTIONS

The following accessories/options are available for your Miami Lumbar system:

- Posterior air bladder for additional lumbar support
- TLSO attachment (includes ATE, PTE, shoulder and axillary straps)
- PTE attachment (includes PTE and axillary straps only)

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