ENGLISH

- MD**
- Medical Device

DESCRIPTION

The device is a silicone pad designed to distribute pressure from socket on prominent bony areas for the residual limb.

INTENDED USE

The device is a prosthetic interface intended to be used as part of a system that replaces a missing lower limb.

Suitability of the device for the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

INDICATIONS FOR USE AND TARGET PATIENT POPULATION

- Lower limb loss, amputation, or deficiency
- No known contraindications

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Caution: The device should not be in direct contact with the skin at any point. Direct skin contact may lead to skin irritation.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Laminated Sockets

- The plaster model has to be dry or isolated.
- PVA bag is pulled over the cast and then left to dry. Remove the backing from the pads and press them in place.
- Pull over the laminating material. Then pull outer PVA bag over and laminate the socket as usual.
- When the plaster model is removed from the socket, the pads will be sticky.
- Remove the adhesive by soaking the pads quite generously with Thinner or rubbing with a cloth soaked in Thinner.
- While still moist, carefully lift off the glue coated polyurethane film from the top of the pads, leaving the fabric covered silicone behind

Thermoplastic Sockets

- The plaster model has to be dry or isolated.
- Remove the backing from the pads and press them in place.
- Thermoform the socket as usual.
- When the plaster model is removed from the socket, the pads will be loose and sticky.
- Remove the adhesive by soaking the pads quite generously with Thinner or rubbing with a cloth soaked in Thinner.
- While still moist, carefully lift off the glue coated polyurethane film from the top of the pads, leaving the fabric covered silicone behind.
- Glue the pads in place inside the thermoplastic socket.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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