

Technical Manual

---

RESOLVE<sup>®</sup> HALO VEST

		3
EN	Instructions for Use	7
DE	Gebrauchsanweisung	12
FR	Notice d'utilisation	18
ES	Instrucciones para el uso	24
IT	Istruzioni per l'uso	30
NO	Bruksanvisning	36
DA	Brugsanvisning	41
SV	Bruksanvisning	46
EL	Οδηγίες Χρήσης	51
FI	Käyttöohjeet	57
NL	Gebruiksaanwijzing	63
PT	Instruções de Utilização	69
JA	取扱説明書	75
ZH	中文说明书	80
KO	사용 설명서	85

2





Chest Circumference at Xyphoid	Posterior Vest Size
60-121 cm (24-48 in)*	ReSolve® Halo Vest - Short
66-137 cm (27-54 in)**	ReSolve® Halo Vest - Tall







\* Waist stabilizers require trimming if measurement is less than 89 cm (35 in)







\*\* Waist stabilizers require trimming if measurement is less than 86 cm (34 in)







3

Anterior Vest Height Adjustment	Length of Anterior Vest (mid-line measurement from the superior edge to the inferior edge)
Bottom Hole (#1)	19.5 cm (7.68 in)
Middle Hole (#2)	21.5 cm (8.46 in)
Top Hole (#3)	23.0 cm (9.13 in)

# ENGLISH

---

Instructions do not supercede physician's orders or hospital protocol. The order at which you put the vest and ring on is left to the physician.

**Caution: Federal law restricts this device to sale by, or on the order of, a physician.**

A complete ReSolve® Halo system includes a ReSolve® Halo Vest used in combination with a ReSolve® Halo Ring or Össur Halo Ring (open or closed back versions).

## INDICATIONS FOR USE

The ReSolve Halo is recommended for use in traumatic or chronic disorders requiring traction along the cervical spine to align the vertebral structures (reduce deformities), maintain reduction, or provide stabilization.

Indications include cervical spine trauma or surgery.

## CONTRAINDICATIONS/ WARNINGS

The ReSolve Halo should not be used under any conditions where pressure to the skull would not be tolerated. When used as part of a traction system (i. e. with weights), the halo should also not be used when the stability of the cervical spine will not permit traction loads.

**Caution:** Emergency access to chest procedures (located on product label on anterior crossplate) should not be used for the purpose of daily hygiene. Repeated flexure of vest may lead to material failure.

This product is for single patient use only.

## MRI CONDITIONAL/WARNING

The Halo is rated MRI Conditional and is safe to use at either 1.5 or 3 Tesla, depending on which type of skull pins are used. Usage at higher magnetic field strength may result in burns to the skin at the pin sites. Consult IFU contained with Halo Ring and Pins for specific MRI safety information.

## PRE-APPLICATION

1. Lay the patient supine.
2. Select posterior vest size by measuring the patient's chest circumference at the xyphoid process (Exception: with pendulous abdomen measure around the widest part of the costal margin) (**Figure 1**). If patient does not fit into these sizes or adjustments are needed, please contact Össur Customer Service for custom vest options. In general patients shorter than 1.75 meters ( 5'9") will be fitted with the short size and patients taller than 1.75 meters (5'9") will be fitted with the tall size. For a properly fitted vest, the chest circumference at the xyphoid process measurement must be used. See sizing chart (**Figure 2**) to choose correct size vest. **Note:** If any trimming of the waist stabilizers or liner is necessary, it is best to complete this prior to fitting the vest. Sharp edges as a result of trimming should be rounded.
3. Maintaining supine position, measure the patient's sternal length to determine anterior vest size. See sizing chart (**Figure 3**) to determine correct anterior vest height setting. Detach the liner from the bottom anterior shell. Adjust height by pushing on the button and sliding the lower anterior vest shell (**Figure 4**).

**CAUTION:** YOU MUST APPLY AND TIGHTEN THE THREADED KNOB TO SECURE THE HEIGHT OF THE ANTERIOR VEST SHELL. FAILURE TO DO SO MAY RESULT IN BODILY INJURY TO THE PATIENT (**Figure 5**).

Trim excess liner if necessary.

4. Lay the anterior vest on patient. The bottom edge should not extend below the costal margin (lowest rib) and upper edge (center) should be below the sternal notch and centered on the chest. Re-adjust anterior vest height if necessary always remembering to secure the height with threaded knob. Remove anterior vest and set aside.

## VEST APPLICATION

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*Application may require one or two assistants. Ensure upright rods are attached to the posterior vest and slightly loosened.*

1. Log roll patient onto their side (**Option 1**) or raise patient 45 degrees from the waist (**Option 2**) maintaining proper c-spine precautions (**Figure 6a and 6b**).
2. Place posterior vest on patient and return patient to the supine position (**Figure 7a and 7b**).

#### **Proper Posterior Vest Fit (Figure 8)**

- A. Upper edge of vest at the level of C7
- B. Ends of crossplate in firm contact with trapezius muscles not on the spine of the scapula
- C. Superior edges of waist stabilizers at level of xyphoid process
- D. Central relief opening in line with spinous processes

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Ensure upright rods are attached to the anterior vest and slightly loosened.

1. Position the anterior vest on patient's chest so that the lateral superior edge is just below the clavicles (**Figure 9**).
2. Wrap waist stabilizers and waist straps around patient and insert them through the loops on the anterior vest (**Figure 10**). Some trimming of the waist stabilizers may be necessary for smaller patients (**Figure 11**). (Be careful not to trim waist stabilizers too short. This can make vest re-application difficult after liner changes.) Sharp edges as a result of trimming should be rounded. Ensure waist straps are in the alignment tabs to prevent straps from slipping off of the waist stabilizers. Tighten waist straps and fasten with hook piece (**Figure 12**).
3. Thread shoulder straps through shoulder loops on the front of the vest. Tighten and secure making sure the plastic shoulder stabilizers are underneath the straps (**Figure 13**). Trimming of the shoulder stabilizers may be necessary depending on the patient size. Alternatively, the plastic shoulder stabilizers may be removed.

**Note: All straps should be tightly fastened without causing patient excessive discomfort.**

4. Secure waist stabilizers by lining up corresponding holes and applying threaded knob (**Figure 14**). **Threaded knob must be applied and tightened to properly secure the waist stabilizers.** For larger patients where the waist stabilizers do not overlap in the front, two threaded knobs will be necessary. It may be helpful to note the number next the hole that the threaded knob is inserted into. This number can be recorded in the patient's chart for future reference (i. e. during liner changes).

#### **Proper Anterior Vest Fit (Figure 15)**

- A. Shoulder pieces in close contact with patient's shoulder
- B. Shoulder straps secure
- C. Vest below sternal notch and centered
- D. Superior lateral vest extensions in firm contact with delto-pectoral groove, no contact on clavicles
- E. Waist stabilizers cover lower costal margins
- F. Threaded knobs are secure
- G. Waist straps secure and in alignment tabs on waist stabilizers

## SUPERSTRUCTURE

Continue to maintain c-spine alignment throughout superstructure application.

1. Align and attach the ring attachment disk and distraction assembly to halo ring (**Figure 16**). Make sure the post on the attachment disk is placed in one of the holes in the "V" portion of the ring.
2. Loosen all upright-transverse clamps and all vest-upright connectors (**Figure 17**).
3. Insert posterior uprights through posterior upright-transverse clamps and tighten. Insert anterior uprights through anterior upright-transverse clamps and tighten (**Figure 18**). Tighten distraction assembly clamp on the transverse rod (**Figure 19**).

4. Lean on anterior vest while tightening anterior vest-upright connectors (**Figure 20**). Then log roll enough to tighten the posterior vest-upright connectors at the attachment points. If necessary, re-tighten shoulder straps, waist stabilizers, and/or waist straps to obtain firm support of vest against sternum.

### **PROPER SUPERSTRUCTURE FIT (Figure 21)**

- A. Transverse bars parallel to ring and in the same plane (alignment may vary in certain conditions).
- B. Anterior uprights rods out of peripheral vision
- C. All bolts tightened (upright-transverse clamps, distraction assembly clamps, vest-upright connectors)

### **DISTRACTION/COMPRESSION**

The distraction assembly is adjustable in the vertical plane to allow for fine-tuning of distraction or compression.

Use the T-handle allen wrench to rotate adjustable height screw located in the distraction assembly (**Figure 22**).

**The adjustable upright is for fine-tuning only.** For gross distraction, loosen all four transverse- upright clamps and distract on the vertical uprights.

### **CAPITAL FLEXION/EXTENSION**

Loosen the ring attachment disk and position the patient as needed. Re-tighten both connections (**Figure 23**).

### **FULL FLEXION/EXTENSION**

Loosen all four upright-transverse clamps, hold patient by transverse rods, and manipulate patient into desired full flexion or extension of the spine. Tighten all four upright-transverse clamps (**Figure 24**).

### **A/P TRANSLATION**

Loosen both distraction assembly clamps and slide the patient anterior or posterior to achieve desired A/P translation. Re-tighten clamps. This technique may be necessary for a patient with kyphosis (**Figure 25**).

## **FINAL CHECKLIST (Figures 8, 15, 21)**

- Bolts tight
- Superior lateral vest Uprights out of peripheral vision
- Transverse bars parallel to ring and all four uprights at the same height (if possible)
- Shoulder pieces in close contact with patient's shoulder
- Shoulder straps secure
- Upper edge of vest within 2 inches of sternal notch and centered on chest
- Anterior extensions in firm contact with delto-pectoral groove, no contact on clavicles
- Lateral borders cover lower costal margins
- Two threaded knobs secure
- Waist straps secure and in alignment tabs on waist stabilizers
- Upper edge of vest midline with C7
- Ends of crossplate in firm contact with trapezius muscles not on the spine of the scapula
- Superior edges of waist stabilizers at level of xyphoid process
- Central relief openings in line with spinous processes

## **OPTIONAL/HELPFUL HINTS**

**Stockinette:** To wick moisture away from the skin, it may be useful to use a body sock or stockinette under the vest. To apply, roll it up from patient's waist to under their arms prior to vest application. Ensure there are no wrinkles or bunching of the stockinette. Slits may be cut in stockinette to feed shoulder and waist stabilizers and prevent stockinette from rolling down during wear. Stockinette should be replaced during liner changes (**Figure 26**).

**Clearing ear space:** Reversing the adjustable upright can bring the transverse rod higher, and clear the space in front of the ear; however there may be height restrictions in the length of the upright rods of the superstructure (**Figure 27**).

## MODIFICATION OF VEST FRAME

Any shape modification to the vest frame should only be done after heating the plastic as follows:

Heat the plastic with a heat gun on both sides of the frame in the desired region evenly until a surface temperature between 125–200°F (approximately 52–93°C) is reached. The heated region of plastic can then be contoured as needed by gently positioning the plastic into the desired shape, and holding this position until thoroughly cooled.

Össur recommends vest frame shape modification in the following instance:

**Kyphotic Curve:** The middle section of the posterior vest can be modified as necessary by heating as described above and then holding in the desired position while the plastic cools. Reversing the direction of the posterior upright rods can accommodate for a kyphotic curve (**Figure 28**). The anterior vest should be set to the shortest position.

## AFTERCARE

**Note:** these instructions are directed towards practitioners/physicians applying a halo, and are not intended as patient care instructions. Refer to Össur Patient Information Manual for patient home hygiene and care instructions.

1. Provide and explain information to patient for daily care of halo.
2. Schedule a 24–36 hour follow up visit post-application to check pin tightness and for signs of complications
3. 3 week visits for liner changes and pin-site inspections should be scheduled. During pin-site inspections, check for redness, discharge, or other signs of infection.

### PROTOCOL FOR POSTERIOR LINER CHANGES

*Per physician's prescription and established protocol only.*

1. Remove the threaded knob from the front of the vest securing waist stabilizers.
2. Place patient prone with chest bolstered slightly to raise face from table/bed.
3. Detach one waist strap and one shoulder strap from the strap loops on the anterior vest, (both on the same side) by loosening. Loosen and detach both posterior upright-vest connectors. Do not remove the uprights from the transverse rod.
4. Hinge open the posterior vest and remove dirty posterior liner. Replace with a clean liner, matching the hook and loop sections (**Figure 29**).
5. Clean skin with hypoallergenic soap and water.
6. Replace and re-attach completely the posterior vest. Hint: It is best to loop the waist strap through the anterior vest before sliding the waist stabilizer back into place. Also, attaching the hook piece at the end of the strap and cutting both on an angle make re-applying the strap easier. Ensure to tighten all straps. Re-attach upright rods to vest and tighten hex bolts on upright-vest connectors.
7. the threaded knob in the front of the vest must be re-inserted to secure the waist stabilizers (**Figure 14**).

To verify proper posterior vest fit, see Proper Posterior Vest Fit checklist (**Figure 8**).

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*Per physician's prescription and established protocol only.*

1. Place patient supine facing upwards.
2. Remove threaded knob securing waist stabilizers. Loosen and detach waist straps, waist stabilizers, and anterior upright-vest connectors. Do not remove the uprights from the transverse rod (**Figure 30**).
3. Remove dirty anterior liner and replace with clean liner, matching the hook and loop sections.
4. Clean skin with hypoallergenic soap and water.
5. Replace and re-attach completely the anterior vest. Ensure to tighten all straps. the threaded knob in the front of the vest must be re-inserted to secure the waist stabilizers (**Figure 14**). Re-attach upright rods to anterior vest and tighten hex bolts on upright-vest connectors.


To verify proper anterior vest fit, see Proper Anterior Vest Fit checklist (**Figure 15**).

## FOR PERFORMING EMERGENCY CPR

Instructions are provided on the anterior portion of the vest. For additional information see written instructions below:

- ① Unscrew bottom knob

- ② Detach waist straps

- ③ Place hand on upper portion of vest and bend lower section towards chest

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**EN – Caution:** Ossur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards is achieved only when Ossur products are used with other recommended Ossur components. If an unusual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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