












Technical Manual

REBOUND® ACL

	3
EN Instructions for Use	5
DE Gebrauchsanweisung	9
FR Notice d'utilisation	14
ES Instrucciones para el uso	19
IT Istruzioni per l'uso	24
NO Bruksanvisning	29
DA Brugsanvisning	33
SV Bruksanvisning	38
EL Οδηγίες Χρήσης	42
FI Käyttöohjeet	47
NL Gebruiksaanwijzing	52
PT Instruções de Utilização	57
PL Instrukcja użytkowania	62
CS Návod k použití	67
TR Kullanım Talimatları	72
RU Инструкция по использованию	76
JA 取扱説明書	81
ZH 中文说明书	85
KO 사용 설명서	89

ENGLISH

LIST OF PARTS

Figure 1

- a. Proximal Thigh Strap
- b. Proximal Thigh Pad
- c. Distal Thigh Liner Strap
- d. Cable Shroud
- e. Cable Covers
- f. Hinge
- g. Proximal Anterior Shin Strap
- h. Shin Pad
- i. Distal Anterior Shin Strap
- j. Upper Frame
- k. Dynamic Tension System (DTS) with Spool
- l. Condyle Pad
- m. Posterior Calf Strap
- n. Cable
- o. Lower Frame

Figure 9

- p. Shear Knob

Figure 10

- q. Dynamic Tension System (DTS) Sticker

Figure 11

- r. Release Mechanism

INDICATION FOR USE

Conditions benefiting from an agonistic force that varies throughout flexion in proportion to ACL tension applied to the thigh and calf in the sagittal plane. These may include:

- Anterior cruciate ligament (ACL) injury
- Reconstruction
- Non-surgical treatment.

CONTRAINDICATIONS

None known.

WARNINGS

- The brace must be prescribed and fit by a qualified healthcare provider.
- Rebound ACL can be ordered Custom Made and is a field adaptable knee brace and should be modified to fit individual patient contour by an orthotist or qualified practitioner.
- This brace is intended for post-surgical or post-injury rehabilitation. It is not intended to be used in moderate or high activity and/or impact sports. A qualified healthcare provider must approve the types of activities that the patient may engage in while using the brace.
- This knee brace is a supportive brace only and is not intended or guaranteed to prevent injury; treatment results may vary.
- This product has been designed and tested for single patient use and is not recommended for multiple patient use.
- Cable tension should only be adjusted with the patient's leg in full extension. Adjusting cable tension in flexion will increase the risk of overloading the knee.
- The patient must not adjust the cable tension. Any adjustments to the

cable tension must be made by the healthcare provider. The brace has been designed to be adjusted to the proper tension using the provided shear knobs. Adjustment using any other means can result in excessive tension and injury to the knee.

- In the event of accidental adjustment to the cable tension, the patient should immediately discontinue use of the brace and contact the healthcare provider.
- To avoid discomfort and potential impairment of circulation, the patient should avoid sitting on the frame for an extended period of time.
- If the device is damaged, or if the patient experiences increased pain, swelling, sensation changes, excessive pressure, skin irritation or any adverse reactions while using this brace, immediately discontinue use of the brace and contact the healthcare provider.

ENSURE THE END-USER AGREES TO THE FOLLOWING

1. Össur accepts no responsibility for any injury which might occur while using Össur knee bracing products.
2. With use of this type of device, it is possible to receive various types of skin irritation, including rashes and pressure points. If the brace causes discomfort or any unusual reactions, discontinue use and contact your health care provider immediately.
3. If you have any questions as to how this knee brace applies to your own specific knee condition, consult your health care provider.
4. This brace is not intended to be used in moderate or high activity and/or impact sports. If you do not agree to the above conditions, contact your brace provider before using this product.

INITIAL BRACE APPLICATION

1. Unfasten all straps (**Figure 2**). Note: Brace comes installed with a 0° extension stop and a 90° flexion stop. Refer to the Range of Motion (ROM) adjustment section if change in ROM is required. Brace comes installed with standard thickness condyle pads. If change in condyle pad thickness is required, refer to Condyle Pad Adjustment section.
2. With the leg in full extension, place the brace on the leg and place hinge center slightly above (approximately ¾") the joint line so the brace can properly seat itself. Note: Ensure that the cable covers are on the outside of the brace and floating free of obstruction. Check both sides of the brace.
3. Position the shin pad over the shin. (**Figure 3**).
4. While holding frame in place, fasten the proximal anterior shin strap by inserting it through the d-ring, then fasten the distal anterior shin strap by inserting it through the d-ring. (**Figure 4**).
5. Next, snugly fasten the proximal thigh strap (**Figure 5**).
6. Detach the hook and then cut anterior side of the distal thigh liner strap to size using cut lines. (**Figure 6**).
7. Reattach the hook tab and then snugly fasten the distal thigh liner strap. (**Figure 7**).
8. Fasten the posterior calf strap through the posterior slot on the upright to affix the hinge in its centered position at midline A/P (**Figure 8**). Note: If the brace is too anterior, loosen the upper shin strap and tighten the calf. If the frame is too posterior, loosen the calf and tighten the upper shin strap.

FIT ADJUSTMENT

Adjusting Pad and strap length:

All straps can be trimmed or adjusted to length. Remember, straps should not restrict circulation or cause cramping during use. For optimal fit and comfort, ensure Strap Pads are centered on the leg when straps are tightened. If a Strap Pad is too long, remove the pad from strap and trim excess with scissors. To adjust the length of a strap, detach the alligator tab at the end of the strap. Adjust the length as necessary and reattach the alligator tab to the strap.

SETTING INITIAL TENSION

1. Two different shear knobs are included:
 - White (standard load)
 - Black (low load)

Note: The standard load (white) is based on a review by Smith et al KSSTA 2013
2. With the patient's leg in full extension, place the desired shear knob in the hole on the spool. Slowly turn the shear knob in the clockwise direction until the shear knob breaks (**Figure 9**).
3. Discard the two pieces of the shear knob.
4. Place the DTS sticker on the spool to cover the hole. This is to prevent the patient and/or anyone other than the healthcare provider from adjusting the spool. (**Figure 10**).
5. Provide Instructions for Use to the patient. Instruct the patient never to remove the DTS sticker and never to adjust the spool.

RE-TENSIONING DYNAMIC TENSIONING SYSTEM

1. Remove brace from leg if patient is still wearing the brace.
2. Unfasten the DTS system from the distal thigh liner.
3. Flip the DTS system so the underside is exposed.
4. Insert a pen into the hole in the mechanism system and apply outward force in the direction of the arrow to disengage the mechanism system from the spool. While holding the mechanism system in the disengaged position, pull the cable out of the DTS, (this will unwind the cable and cause the spool to spin in the counter clockwise direction) until the DTS slides freely on the cable (**Figure 11**).
5. Center the DTS by pulling the DTS away from the brace and centering it between the uprights, (ensure that the cable is equal length on each side of the DTS) (**Figure 12**).
6. Place the desired shear knob into the hole on the spool and slowly turn the shear knob in the clockwise direction approximately a half turn to lock the DTS onto cable.
7. Flip the DTS back over and re-affix in its original position on the distal thigh liner. Note: The cables should be loose but of equal length on each side of the DTS.
8. Repeat steps 1–8 in the "Initial Brace Application" instructions and steps 1–5 in the "Setting Initial Tension" to re-tension the DTS.

RANGE OF MOTION ADJUSTMENT (ROM)

Adjustment should be performed by a health care professional. The brace comes with 0° extension stops and 90° flexion stops attached to the hinges. If needed, ROM control can be adjusted using the extension and flexion stops.

1. Remove the condyle pads from the medial and lateral hinges. Remove retaining screws. Using a screwdriver or pen, slide the 90° stop out of the hinge on both sides.
2. Re-apply the desired flexion stop and gently secure it with the retaining screw (**Do not over tighten**).
3. Follow the same procedure to change the extension stops.

CONDYLE PAD ADJUSTMENT

Attached to the hinges are the standard size pads. Pairs of thicker and thinner pads are included. The pads should lightly touch the medial and lateral condyles of the knee.

Widening or narrowing fit at knee joint: the condyle pads should lightly touch the medial and lateral condyles of the knee. To widen or narrow the width at the knee, the medium condyle pads can be removed and replaced with the thicker or thinner version shipped with the Rebound ACL. Any combination of pads can be used to achieve a more comfortable fit.

BRACE MAINTENANCE

The brace is designed for minimum maintenance and care. Follow these easy maintenance steps to keep the brace in good condition.

CLEANING

- Liners, straps and pads may be removed for cleaning
- Do not machine wash
- Hand-wash in luke warm water (30°C) using mild non-detergent soap
- Rinse thoroughly
- Do not place in dryer. Only air-dry. If migration occurs, clean silicone liners with an alcohol wipe or damp cloth. If brace is used in salt water or chlorinated water, rinse well and air dry.

HINGE MAINTENANCE

- Do not disassemble the hinge
- Check regularly for foreign materials in hinge (i. e. dirt or grass)
- Carefully clean dirt from hinge and cable covers using water or compressed air
- Lubricate hinge with lithium grease

ACCESSORIES

Accessories are available for your Rebound® ACL brace at an additional charge. To order call your Össur Customer Service Department.

DISPOSAL

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations. Users should contact their local governmental office for information on how these items can be disposed of in an environmentally sound manner.

PRODUCT MATERIAL LIST

Aluminum, Acetal, ABS, Brass, EVA, Glass-filled Nylon, Nylon, Polyurethane, Polyester, Stainless Steel, PTFE.

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