

Instructions for Use

---


VARI-FLEX® MODULAR


		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	12
FR	Notice d'utilisation	19
ES	Instrucciones para el uso	26
IT	Istruzioni per l'uso	33
NO	Bruksanvisning	40
DA	Brugsanvisning	46
SV	Bruksanvisning	52
EL	Οδηγίες Χρήσης	58
FI	Käyttöohjeet	65
NL	Gebruiksaanwijzing	71
PT	Instruções de Utilização	78
PL	Instrukcja użytkowania	85
TR	Kullanım Talimatları	92
RU	Инструкция по использованию	98
JA	取扱説明書	105
ZH	中文说明书	111
KO	사용 설명서	117






2

3

### DESCRIPTION

The device is a prosthetic foot. It is composed of a base foot, with carbon-fiber heel and keel.

The device consists of the following components (**Fig. 1**):

- A. Foot Module
- B. Heel
- C. T-Adapter
- D. Heel Wedge
- E. Male Pyramid

This device must be used with a Foot Cover and a Spectra Sock.

### INTENDED USE

The device is intended as a part of a prosthetic system that replaces the foot and ankle function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

### *Indications For Use and Target Patient Population*

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low to high impact use, e.g., walking and occasional running.

### GENERAL SAFETY INSTRUCTIONS

**Warning:** Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

### DEVICE SELECTION

Verify that selected variant of the device is suitable for the impact level and weight limit according to the following table.

**Warning:** Do not exceed weight limit. Risk of device failure.

Incorrect category selection may also result in poor device function.

Weight kg	117-130	131-147	148-166	167-185	186-203
Weight lbs	257-287	288-324	325-356	366-408	409-448
<b>Low Impact Level</b>					
Size 22-24	7	8	9	10	11
Size 25-27	N/A	N/A	N/A	10	11
Size 28-30	N/A	N/A	N/A	10	11
<b>Moderate Impact Level</b>					
Size 22-24	8	9	10	11	N/A
Size 25-27	N/A	N/A	10	11	12
Size 28-30	N/A	N/A	10	11	12
<b>High Impact Level</b>					
Size 22-24	9	10	11	N/A	N/A
Size 25-27	N/A	10	11	12	N/A
Size 28-30	N/A	10	11	12	N/A

## ASSEMBLY INSTRUCTIONS

### Foot Cover and Spectra Sock (Fig. 2)

**Caution:** To avoid pinching fingers, always use a shoehorn.

1. Put the foot into the Spectra Sock.
2. Remove the foot cover attachment.
3. Use the applicable (straight) end of a shoehorn to put the foot with the Spectra Sock into the Foot Cover.
4. Move the shoehorn up to fully push the foot into the Foot Cover.
5. Put the foot cover attachment back onto the Foot Cover.
6. After alignment is complete, fix the Spectra Sock to the prosthesis to seal against dust and dirt.

**Note:** The Spectra Sock must be pulled up to prevent it from interfering with moving parts of the foot.

If required doff the Foot Cover as follows:

1. Remove the foot cover attachment.
2. Insert the applicable (edged) end of a shoehorn behind the foot.
3. Push the shoehorn down and pull the foot out of the Foot Cover.
4. Fully remove the Spectra Sock.

### Bolt Assembly

To attach a Pyramid or a Tube, thread each bolt through washer and Foot Module into Pyramid/Tube. Apply Loctite 410 to the bolt and torque according to the relevant instructions. When using a Tube, place friction pad between foot and Tube.

### Cutting Foot Module (Fig. 4)

Before cutting check for enough clearance

- Measure for length and add 25 mm for T-Adapter, 38 mm when using a Pyramid.

**Warning:** The T-Adapter/Pyramid should only be placed where it makes full contact with the carbon blade.

- Clamp Foot Module in vise. Use a sharp fine tooth hacksaw blade.
- Do not cut below the cut limit sticker line for attachment.

- Sand away any sharp edges.
- Secure foot to Pyramid or T-Adapter using C-clamp
- Check height of prosthesis on patient while standing between parallel bars.
- Cut to final length.

### ***Drilling the Foot Module Attachment Holes***

- Back-up Foot Module when drilling using wood block or discarded carbon fiber section and clamp securely.
- Use drill guide template to locate holes (**Fig. 5**).
- Use sharp drill bits. Bullet nose or solid carbide tipped drill bits suggested.
- Remove template and drill to appropriate hole size.

Drill size:

Cat 1-6: 8.5 mm

Cat 7-9: 10 mm

- Drill speed 1200-1800 RPM.
- Feed drill slowly with light intermittent pressure - do not force.
- Foot is now ready for assembly.

### ***Using Modular Alignment Adapter***

- Attach the Alignment Adapter to the Foot Module (**Fig. 6**).
- Make sure grit is between Alignment Adapter and Foot Module.
- Connect the foot and the Alignment Adapter to the prosthesis.
- Adjust the Alignment Adapter to the correct height by sliding it on the blade.
- Clamp the Alignment Adapter to the Foot Module by tightening the screws to 22 Nm.
- Disconnect the Foot Module, with the Alignment Adapter clamped on to it, from the prosthesis.
- Drill appropriate holes according to **Drilling Foot Module Attachment Holes** instructions through the drill bushings instead of the guide template. Note that the drill bushings can only accommodate 8.5 mm drills.
- Therefore a second drilling with a 10.5 drill is required for the larger categories.
- Mark or measure correct length of Foot Module.
- Detach Alignment Adapter from the Foot Module.
- Cut Foot Module according to **Cutting Foot Module**.
- Attach Pyramid to Foot Module.
- Connect the Foot Module with the Pyramid attached to the user.

### ***Lamination Connector***

The Lamination Connector must be laminated directly into the posterior wall of the socket.

**Warning:** Before starting make sure to orient the Lamination Connector the correct way. The shoulders of the inserts (marked with arrow) should be oriented towards the socket (**Fig. 7**).

1. Put the Lamination Connector in the middle third of the socket to make height adjustment easier.
2. Make sure the socket surface is flat. A flat surface below the Lamination Connector distributes the load to the connected foot.
3. Use a friction pad between the Lamination Connector and the foot.

## **Prosthesis**

Assemble prosthesis with applicable devices.

**Warning:** Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

**Warning:** Ensure proper attachment by following the applicable device assembly instructions.

Remove the protective film on the pyramid after fitting.

## **ALIGNMENT INSTRUCTIONS**

### **Bench Alignment (Fig. 3)**

#### **Alignment Goal**

Alignment reference line (B) should:

- pass through midpoint of socket at either the patellar tendon level or at the ischial tuberosity level (D).
- fall at the 1/3 mark on the inside of the foot cover.

**Note:** Prioritize knee alignment over foot alignment if there is a mismatch.

#### **Alignment Instructions**

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect either the socket or the knee to the foot and establish the correct knee center height.
3. If using a prosthetic knee: Position knee according to knee alignment instructions
4. On the lateral side of the socket, make a first mark (alignment reference point) at the midpoint of the socket at either patellar tendon level or the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through the alignment reference point (D).
6. Set the correct socket angles for flexion/extension and abduction/adduction.
7. If using a prosthetic knee: Use the applicable adapters to connect the socket to the prosthesis.

#### **Static Alignment**

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

#### **Dynamic Alignment**

Ensure the patient is familiar with the functioning of the device.

The heel to toe action can be influenced by:

- Heel Stiffness
- Anterior-Posterior positioning of device
- Dorsi-Plantarflexion
- Shoe Characteristics

Consider the following actions if needed:

#### **Symptoms**

- Device comes to flat position too early (patient feels they are sinking into a hole)
- Climbing over the toe requires extra energy

- Toe feels too stiff
- Knee hyperextends

#### **Action**

- Check Wedge instructions (if applicable)
- Shift socket anterior (or device posterior)
- Consider dorsiflexion
- Check heel of the shoe and shoe performance

#### **Symptoms**

- Rapid heel to toe movement
- Poor control over prosthesis at initial contact
- Minimal energy return feeling
- Too little push off from the toe
- Knee becomes unstable

#### **Action**

- Check Wedge instructions (if applicable)
- Shift socket posterior (or device anterior)
- Consider plantarflexion
- Check heel of the shoe and shoe performance

#### **Heel Wedges (Fig. 8)**

The small, medium and large Wedges are used to change heel stiffness. The Wedges can be trimmed using sharp scissors to customize stiffness. For temporary Wedge placement, use tape to secure the wedge in position.

#### **For permanent Wedge placement**

- Roughen the upper and lower surface of the heel with abrasive paper.
- Apply instant adhesive on the lower side of the Wedge only.
- Locate in the foot/heel junction and position before adhesive sets.
- If the foot has a split toe, remove a thin slice in the middle by cutting with a sharp knife through the split in the carbon blades.

For removal, the adhesive may be softened by soaking in acetone or cyanoacrylate adhesive remover.

#### **USAGE**

##### ***Cleaning and care***

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

##### ***Environmental Conditions***

The device is Waterproof.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt.

Continuous exposure is not allowed.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

Dry with a cloth after contact with fresh water or humidity.

#### **MAINTENANCE**

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

### Noise from Foot

Noise may occur if sand or debris is present in device. In that case, the healthcare professional should doff the foot, clean it with the help of compressed air and replace the Spectra sock if it is damaged.

### REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

### DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

### LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

### Compliance

This device has been tested according to ISO 10328 standard to two million load cycles.

Depending on patient activity this may correspond to 2-3 years of use.

ISO 10328 - Label		
Category	Weight (Kg) Moderate Impact	Label Text
7	116	ISO 10328-P6-116 kg
8	130	ISO 10328-P7-130 kg
9	147	ISO 10328-P8-147 kg
10	166	ISO 10328-P8-166 kg
11	185	ISO 10328-P8-185 kg
12	203	ISO 10328-P8-203 kg
*Actual test loads reflect maximum body mass		

ISO 10328 - "P" - "m"kg \*)

\*) Body mass limit not to be exceeded!

For specific conditions and limitations of use  
see manufacturer's written instructions on  
intended use!

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