




Technical Manual














i - DIGITS














3

EN		Technical Manual	6
CS		Technický manuál	21
SK		Indikácie použitia	37
RO		Indicații de utilizare	53
HR		Upute za uporabu	70
HU		Használati javallat	86
BG		Указания за употреба	103
SL		Indikácie pre použitie	120

ENGLISH

INSTRUCTIONS FOR USE

i-DIGITS

The i-Digits is referred to as “*the device*” in the following document. This document provides information on the indications for use and handling of the device. It is intended for a certified clinician. The device may only be fitted and configured by a qualified practitioner authorized by Össur after completing the corresponding training. These “*Technical Manual*” relate to: i-Digits Quantum, i-Digits Access, i-Digits France, analog electrodes and mains charger.

PRODUCT DESCRIPTION

i-Digits is a fully customized partial hand prosthesis with individually powered digits. i-Digits work in conjunction with remaining fingers (**Fig. 1a**). The wristband houses the microprocessor and removable batteries. It is fastened using a buckle and is designed to fit loosely around the distal forearm. The flexible covering/ bellows allows access to wiring during prosthesis assembly. The product label can be found on the wristband (**Fig. 1b**). i-Digits Quantum has 20 different grip options available and 12 customizable my grips. All other i-Digits devices have 12 grip options available. The grips can be accessed and programmed through the Biosim or My i-Limb app. The available grips for the device are displayed on the quick grips page of the Biosim or My i-Limb app.

Features comparison			
Control options	i-Digits Quantum	i-Digits France	i-Digits Access
App control	Yes	Yes	Yes
Muscle control	Yes	Yes	Yes
Proximity control	Yes	Yes	-
Gesture control	Yes	-	-
Grips available	20	12	12
my grips	12	-	-

INDICATIONS FOR USE

- Upper limb amputation
- Congenital upper limb absence

CONTRAINDICATIONS

None known.

INTENDED USE

The i-Digits are intended as a part of a prosthetic system that replaces function of a missing upper limb.

REQUIRED DEVICES

The Biosim and the My i-Limb mobile app requires an Apple iOS device supported by the manufacturer, e.g. an iPhone or iPad. See the Apple Store for device compatibility.

SAFETY INSTRUCTIONS

Warnings

i-Digits:

- The end user is the intended operator of the device and is responsible for its use.

- Improper handling or adjustment of the i-Digits may cause malfunction of the device.
- The i-Digits do not provide sensation, heat and moisture cannot be felt. The i-Digits is for low to moderate impact activities only.
- Do not use without approved cover(s).
- Do not use with damaged cover(s).
- Do not disassemble componentry or modify in any way.
- Do not service or perform maintenance when in use.
- Do not carry objects using only the tips of the digits. Carry objects by evenly distributing weight across the digits, close to the knuckles and palm of the hand (**Fig.2**).
- Do not use for heavy lifting.
- Do not use with machinery with moving parts that may cause personal injury or damage.
- Do not use for extreme activities that may cause injury to a natural hand.
- Do not expose to vibration.
- Do not expose to excessive or high forces, particularly on the fingertips and on the side of the digits.
- Do not expose to water.
- Do not expose to excessive moisture, liquids, dust, high temperatures, or shock.
- Do not use in hazardous environments.
- Do not expose to flames.
- Do not use in or expose to explosive atmospheres.
- The electrode is an APPLIED PART.
- This device is for single patient use.
- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.
- The electrodes may contain nickel.

Batteries:

- Do not bend or exert excessive pressure on the battery.
- Do not pierce the battery.
- Do not disassemble battery.
- Do not expose batteries to high temperatures.
- Do not incinerate batteries.
- Do not short circuit the battery.
- Do not store batteries inside a vehicle.
- Dispose of batteries in accordance with US, European or local regulations.

Precautions

i-Digits:

- Users must comply with local regulations on the operation of automobiles, aircraft, sailing vessels of any kind and any other motorized vehicle or device. It is entirely the user's responsibility to seek confirmation that they are physically and legally able to drive using the i-Digits and to the fullest extent permitted by law.
- Only use with approved Össur accessories and tooling.
- Maintenance, repairs and upgrades may only be performed by qualified Össur technicians and technical partners. Össur will provide upon request information to assist service personnel in repair of device.
- Using an i-Digits device to operate electronic devices connected to a mains outlet may impact function.
- It is not recommended to use the device while in close proximity of other medical electrical equipment.
- Do not use while charging is in progress.
- Only use with covers supplied by Össur.
- Always use with covers to avoid risk of electrostatic build-up and discharge.
- Do not use oil-based lotions on the skin, e.g. Vaseline.
- Do not expose electrode to dirt or fluids.

Batteries:

- Only use Össur batteries with this device.
- Only use the Össur charger to charge Össur batteries.
- Make sure the battery is not subject to continued pressure once fitted.
- Batteries are to be replaced annually.

- If the battery has visibly ballooned or swelled:
 - discontinue the charging process immediately
 - disconnect the battery
 - remove to a safe area
 - leave and observe for 15 minutes
 - replace with new battery
 - do not re-use
 - dispose of any leaking batteries in an appropriate manner

Use of the Battery charger should cease if any of the following occur:

- Battery swelling/explosion/fire.
- Charger power supply burnout/fire.
- Red LEDs at each battery cradle indicating battery fault.
- Batteries, charger base, or power supply surfaces become hot.

Do not position the mains power supply to make it difficult to disconnect the device from the mains power.

If the device or charger will not be used for a long period of time, it is advisable to remove the battery from the prosthesis.

POWER

Batteries

The device can be employed for up to 600 full opening and closing cycles from a single charge of a set of 800mAh batteries. The number of open and close cycles achieved will vary depending on battery age and device use. After one set is drained, the batteries can be easily swapped to the back-up set.

Batteries are placed into the battery case on the wristband, where they are secured by the holding catch located towards the wrist. Release the holding catch to remove the battery (**Fig. 3**).

Each battery is equipped with an LED indicating the charge level of the battery. When the battery charge level is low, the LED will light up red. The LED will stay lit until the battery has been charged to an adequate level.

Charging the Batteries

Charging time: 90 minutes to 3 hours.

Remove the batteries from the prosthesis and insert into the charger base unit. Connect the charger base unit with the power cable. Plug the power cable into the power outlet.

The charging state is indicated on the back of the base unit (**Fig. 4**):

- Middle light on: charger is plugged in
- 2nd and 5th lights flashing green: batteries are charging
- 2nd and 5th lights solid green: batteries are charged
- 1st and 4th red lights on: battery fault, unplug and try again. If lights continue to illuminate, the user should contact the clinician.

Turning On/ Off

Turn the device on by pressing the wristband button; the LED will light up for several seconds to indicate power is on. To turn off, press the button again; the LED will light up for a short moment.

CONTROLLING THE I-DIGITS

There are various control options for accessing automated grips on i-Digits devices. Control options vary between different models.

Identifying the Device Number:

Each i-Digits has a unique device number for identifying the device located on the wristband underneath the battery (**Fig. 5**).

When connecting the i-Digits device to the Biosim or the My i-Limb app, the device number will be displayed on

the connection screen. Selecting the number connects the app to the i-Digits device. Alternatively, when connected to the app the device number is displayed in the 'about' section.

Gesture Control (available only on i-Digits™ Quantum)

Enables an automated grip to be accessed through a smooth motion of the prosthesis in one of four directions (forward, back, left or right). The grips programmed for each direction are customized to the user's requirements using the Biosim or My i-Limb app.

To access gesture control:

- Device must be in full open hand mode.
- Hold arm parallel to the ground (elbow bent to 90°).
- Maintain an open signal until finger twitches.
- Move the hand within 1 second in direction set to access a desired grip.
- i-Digits Quantum will adopt the grip.

NOTE: Gesture control can also be accessed using co-contraction, depending on the set-up of the i-Digits through the prosthetist.

App Control

An automated grip can be accessed by touching an icon in the My i-Limb app. These automated grips are called quick grips. i-Digits will exit the grip when the icon is tapped again or by selecting another grip icon to enter.

Muscle Control

Triggers are specific muscle signals that can be used to access an automated grip. Four triggers are available: hold open, double impulse, triple impulse and co-contraction.

Muscle control can be activated and programmed using the app.

Proximity Control (not available on i-Digits Access)

Grip Chips are small Bluetooth devices which activate grips on the prosthesis when the prosthesis is moved close to them or by tapping the Grip Chip.

Before operating a Grip Chip for the first time, pull the plastic tab out from the device to engage the battery.

To use Grip Chips, make sure the i-Digits device is not connected to the Biosim or My i-Limb app.

- Proximity is accessed by fully opening the i-Digits near (15 cm/ 6") to a Grip Chip and then relaxing the muscles. Wait until the grip is accessed, which may take up to 3 seconds. To exit a grip accessed via proximity control, give a long open signal.
- Double tapping a Grip Chip activates the grip. Quickly tap a Grip Chip twice, just like double-clicking a computer mouse. The LED on the Grip Chip flashes once when tapping is successful.

NOTE: A 3 second pause between each double tap is required. The pause prevents the grip chip from incorrectly detecting multiple taps in a very short period of time. This can cause the hand to enter a grip and then immediately exit it again.

The i-Digits must be fully open/digits stalled open for proximity or tapping to be successful.

The individual Grip Chips are programmed using the Biosim and the My i-Limb app and can be reprogrammed by the user at any time.

Replacing the Grip Chip battery

The battery is expected to last 6 months. Battery life is reduced if Grip Chip is stored close to the i Digits device. Battery is not rechargeable.

- To replace the battery, insert fingernail between upper and lower caps and run around joint to release.
- Remove battery. Insert new battery under retaining clip, ensuring writing on battery is facing outwards.
- Replace the caps by pressing upper and lower caps together.

Battery type: CR1616. Replacement batteries are available from customer support upon request.

Wristband, Dummy Positioning and Digit Assembly

Positioning the Socket Terminal Dummy

Place the socket terminal dummy on the dorsal surface of the positive cast. The position should be central to the midline of the wrist joint. Leave 30 mm space between the socket terminal dummy and the center of the wrist joint, to enable good range of motion at the user's wrist (**Fig. 6**).

Digit Removal from the Knuckle

Use the T8 screwdriver to remove the M3 grub screw from the knuckle assembly (**Fig. 7**). The digit can now be freely disengaged from the knuckle assembly.

Extra Small and Medium Digit Assembly

Depending on the size of the kit, the knuckle mounting plate will vary in shape.

Figure 8 shows a medium size digits assembly.

Figure 9 shows an extra small size digits assembly.

Knuckle Mounting Plate

The knuckle mounting plate supports 4 digits. Cut to adjust to the number of digits required by the user (**Fig. 10**).

Attaching the Knuckle Mounting Plate

Use the T6 screwdriver to attach the knuckle assembly to a knuckle mounting plate using three M2x6 mm Torx screws (**Fig. 11**).

Lamination tabs can be bent or shortened as appropriate.

Consider fabrication before cutting and removing tabs, since some tab length will be needed to fix the plate to the frame.

Digit Attachment to the Knuckle Assembly

Pre-flex the digit slightly to aid assembly.

Firmly place the digit into the knuckle assembly.

Insert M3 grub screw using T8 bit in torque driver. Torque to 1 Nm (**Fig. 12**).

NOTE: Screws are for single use only.

Alignment Transfer

Thumb Alignment

Use the thumb alignment plate to help position the thumb mounting plate (**Fig. 13**).

Screw the thumb alignment plate onto the index position of the knuckle mounting plate using three M2x5 mm screws.

Depending on if it is a left- or right-handed device, attach the thumb rotation assembly to the correct side of the alignment plate using three M2x5 mm screws (**Fig. 14**).

Knuckle Dummies

Use single knuckle dummies to help transfer alignment. Attach the single knuckle dummies to the knuckle mounting plate using three M2x16 mm countersunk Torx screws and the T6 screwdriver (**Fig. 15**).

Insert the alignment bars into the single knuckle dummies. Lock each in place using M2x5 mm countersunk torx screws and the T6 screwdriver (**Fig. 16**).

Control Inputs

Remote Electrode Domes

Assemble the electrode domes into the silicone socket by integrating the suspension piece into the silicone. Use the dummies supplied to create optimal spacing and relief (**Fig. 17**).

Check Socket Electrodes

Use the check socket electrodes during the check socket process (**Fig. 18**). This allows the position of the electrode to be moved during the test fitting stage.

Compact Electrodes

During fabrication of the silicone socket, use the supplied dummies to make a recess at the correct position for the compact electrode. The contact surface of the compact electrode must protrude beyond the inner surface of the silicone socket to ensure skin contact.

Force Sensing Resistors (FSRs)

When using FSRs, the FSR (**Fig. 19**) should be positioned between the silicone socket and the inner lamination. A flat surface must be created on the inner lamination to ensure the FSR is not bent when assembled.

Finishing Options

Two finishing options are possible, with knuckle fairings either incorporate or exposed within the lamination. If using incorporated fairings, the supplied modified fairings and knuckle dummy must be used during lamination. See **Figs. 20 and 21**.

Assembly

Thumb Rotator Assembly – Setting and Assembly

Assemble the thumb rotator in the correct order of parts (**Fig. 22**). To increase the force required to rotate the thumb use the 8 mm spanner to tighten the M4 locking nut. The grub screw will need to be loosened to tighten the nut.

NOTE: make sure the clutch plate is not contaminated as this will affect performance.

Use fingers to tighten the nut and then turn a further 120° to 150° using the 8 mm spanner to set thumb friction to patient preferences (**Fig 23**). The wire will pass through the center of the spigot.

Attaching the Thumb

Remove the base cap on the thumb rotator assembly by removing the two M2x5 mm screws (**Fig. 24**).

Remove the knuckle fairing cover from the thumb assembly by removing the two M2x4 mm screw. (**Fig. 25**).

Remove the thumb from the knuckle assembly by removing the M3 grub screw (**Fig. 26**).

Attach the knuckle assembly onto the thumb rotator assembly using three M2x5 mm screws (**Fig. 27**).

Put the thumb wire through the center of the thumb pivot.

Use a torque driver set to 1 Nm to attach the thumb to the single knuckle assembly (**Fig. 28**). Make sure that the thumb wire is not caught between the cap and the pivot assembly.

Move the knuckle fairing over the rear of the thumb. Use two M2x5 mm screws to attach the knuckle fairing.

Re-attach the base cap, using the screws previously removed (**Fig. 29**).

Insert the thumb digit wires into the provided connector. Make sure digit wire 'tooth' is fully secured into the connector (**Fig. 30**).

After diagnostic fitting, remove the thumb wire connector to enable fabrication of the definitive prosthesis.

To remove connector from thumb wire:

Place small flat screwdriver under tab of connector and lift up (**Fig. 31**).

Gently pull on digit wire to release from connector.

NOTE: use a new connector each time.

Wristband Wiring Connections

Remove the bellows cover by gently pulling it over the socket end of the bellows assembly.

Remove bellows from wristband (**Fig. 32**).

Remove screws to release top collar and PCB from bellows assembly (**Fig. 33**).

Connecting the Digits to the Wristband

(**Fig. 34**)

Slide the digit wires through the aperture at the front of the bellows assembly.

Insert the digit wire connectors into the bellows PCB. A click should be felt when fully inserted (**Fig. 35**).

Create slack in the digit wires, so that the bellows assembly flexes without causing tension in the digit wires.

(**Fig. 36**)

NOTE: if there is not enough slack, flexing the bellows may lead to failure of the digit wire connection.

Removing the Digit Wires

To remove a digit wire, pull the connector piece to remove the digit wires of a specific digit.

Connecting the Electrodes to the Wristband

When digit wires are connected, slide electrode cables through the aperture at the distal end of the bellows (**Fig. 37**). Connect the electrodes cables as shown in the mapping (**Fig. 34**).

Mate the connectors by aligning the mating axes of both connectors and then pushing them together. An audible “click” will confirm a fully mated connection.

Route the signal cables inserted in signal 1 position to avoid excessive bending of the connector (**Fig. 38**).

Check cables are tension-free when bellows is fully flexed (**Fig. 36**).

Check that digits and electrodes are all functioning correctly before completing the next stage.

Removing the Electrode Cables

Use the extraction tool supplied to remove electrode cables when necessary.

Insert the end portion of the tool under the connector flanges (**Fig. 39**) and pull off vertically, in the direction of the connector mating axis.

FSR PCB

If using an FSR wristband, insert the FSR cables through the aperture at the socket end of the bellows (**Fig. 40**).

Insert the pairs of wires into the FSR receivers using the same method as inserting the digit wires. Ensure click is heard indicating secure connection.

Remove the FSR wires the same way as removing the digit wires.

Completing the Wristband Assembly

Carefully reposition the wired PCB into the upper and lower collar. Screw into position to secure (**Figs. 41**).

NOTE: take care not to damage the PCB.

Slide the bellows cover back into position over the bellows.

Bond the bellows cover to the bellows at either end using a few spots of superglue.

Reattach the bellows by pushing the connector from the bellows into the connector from the wristband (**Fig. 42**).

Insert the two M2.5 screws (**Fig. 43**).

Attach the wristband to the socket using the M4 socket button head screw.

Digit Covers

Attaching the Digit Covers

A digit cover must be used for each digit. Digit covers can be slid onto a digit.

The same applies for the thumb and thumb cover.

Make sure covers are fully to the tip of each digit.

CLEANING

The i-Digits wristband, digits, cover and electrode surface can be cleaned with a soft damp cloth and mild soap.

Clean the electrode surface regularly.

Do not use any strong chemicals.

NOTE: Do not submerge the i-Digits, wristband, cover, electrode surface or charger in water.

Clean the cover with isopropyl alcohol once a week to help with disinfection.

MAINTENANCE

Instruct user to always turn off the i-Digits device when not in use.

Instruct user to charge the battery after use.

A certified clinician can perform limited maintenance and repair activities. Do not perform any maintenance or repairs which are not described in this manual.

Make sure the i-Digits device is returned to Össur for service every 12 months.

TROUBLESHOOTING

Problem	Action
Prosthesis does not operate	Ensure prosthesis is turned on
	Ensure battery is connected
	Ensure battery is charged
	Ensure electrodes have good contact
One digit does not operate	Use Health Check in the Biosim app to check if the digit is operating correctly
Prosthesis stops halfway during an action	Electrode settings may need to be adjusted
	Check the electrode cable is not damaged
	Check the digit screws are not loose
	Ensure electrodes have good contact
User complains the prosthesis is difficult to operate	Electrode settings may be too low, review on Biosim myograph
	Ensure the battery is sufficiently charged
	Ensure the electrodes are grounded properly and maintain contact with the user's skin throughout use
	Check electrode placement and wiring
Battery charge does not last as long as expected	Fully charge the battery overnight. Check the battery connection.
	Ensure electrodes are not set too high
	Check user is not holding a sustained signal to the hand by reviewing real time graph, or check for excessive signal counts in usage statistics
	Replace with a spare battery
Battery is not working	Check the battery is connected
	Check the battery is charged
	Check the device using one of the spare batteries

TECHNICAL SPECIFICATIONS

i-Digits	
Voltage	7.4 V (nominal)
Max. current	5 A
Battery capacity	Rechargeable lithium polymer 7.4 V (nominal); 800 mAh capacity
Max device load	20 kg/ 44lbs
Finger carry load (static limit)	5 kg/ 11lbs
Time from full open to full close	0.8 seconds
Expected service life	5 years
Power Supply Specifications	
Manufacturer	FRIWO
Model no.	FW8030M/24
Input	100–240 VAC, 50–60 Hz, 0.6–0.3 A
Output	24 VDC, 1.25A
Charging Base	
Model no.	PL069653
Input	24VDC, 1.25A

ENVIRONMENTAL CONDITIONS

Do not use, transport or store the i-Digits outside of the boundaries in the table below:

	Use	Shipping	Extended storage
Temperature	0°C to 40°C (32°F to 104°F)	-40°C to 70°C (-40°F to 158°F)	-25°C to 70°C (-13°F to 158°F)
Relative humidity	0% to 95%	0% to 95%	10% to 95%
Atmospheric pressure	700 hPa to 1060 hPa	700 hPa to 1060 hPa	700 hPa to 1060 hPa

COMPLIANCE

Device has been tested by an independent accredited laboratory and found to comply with:

- 1. EN60601-1:2006/A1:2013
- 2. EN60601-1-2:2015

Charging Base:

The device complies with:

- IEC 60601-1 ed. 3.1 Tested at a CB-certified CSA Canada Laboratory.

Related Standards:

- IEC 60601-1, IEC 60601-1-6 & IEC 60601-1-9
- Protection against electrical shock – Class II ME Equipment
- Degree of protection against electrical shock – Type BF provides additional protection against electric shock
- Degree of protection against ingress of water – IP2X (**IEC 60529**)
- EMI/EMC
- Compliance with IEC 60601-1-2:2014 +A1: 2020 + Canadian national deviations, EN 60601-1-2/A1: 2021

BLUETOOTH MODULES REGULATED INFORMATION

This device contains the following radio frequency transmitters:

Model	Re	Type and Frequency Characteristics	Effective Radiated Power
Bluetooth Low Energy Dual Mode Module Model BR-LE4.0-D2A	FCC Contains FCC ID: XDULE40-D2 Canada Contains IC: 8456A-LE4D2 Japan Contains transmitter with certificate number  R 205-160268	(Dual Mode) Version V2.1 +ED (GFSK + $\pi/4$ DQPSK + 8DPSK) 2402 -2480 MHz Version V4.0 (GFSK) 2402 -2480 MHz	Adjustable Power (-23 dBm to 10.5dBm) short to long range

ELECTROMAGNETIC COMPATIBILITY

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity and result in improper operation.

In order to regulate the requirements for Electromagnetic Compatibility (EMC) with the aim to prevent unsafe product situations, the BS EN 60601-1-2/ IEC 60601-1-2 standard has been implemented for all Össur Myoelectric Prosthetic Devices. This standard defines levels of electromagnetic emissions for medical devices.

The device and charging base are suitable for use in any environment except where immersion in water or any other fluid is possible, or where exposure to high electrical and/or magnetic fields can occur (e.g. electrical transformers, high-power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

Refer to further guidance below regarding the EMC environment in which the device should be used:

i-Digits Guidance and manufacturer's declaration – electromagnetic emissions		
Össur Myoelectric Prosthetic Devices are intended for use in the electromagnetic environment specified below. The customer or the user of the Össur Myoelectric Prosthetic Devices should assure that it each are used in such an environment		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1 – Class B	Össur Myoelectric Prosthetic Devices use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. Össur Myoelectric Prosthetic Devices are suitable for use in all establishments, including domestic establishments and those directly connected to the public low- voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Group 1 – Class B	
Harmonic emissions IEC 61000-3-2		
i-Digits	Not applicable - Battery powered	
Charging Base	Group 1 – Class B	
Voltage fluctuations/ flicker emissions IEC 61000-3-3		
i-Digits	Not applicable - Battery powered	
Charging Base	Group 1 – Class B	

Guidance and manufacturer's declaration – electromagnetic immunity			
Össur Myoelectric Prosthetic Devices are intended for use in the electromagnetic environment specified below. The customer or the user of the Össur Myoelectric Prosthetic Devices should assure that each are used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV power supply lines ± 1 kV input/output lines	i-Digits: Not applicable Charging base: ± 2 kV power supply lines ± 1 kV input/output lines	Not applicable Battery Powered
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	i-Digits: Not applicable Charging base: ± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Not applicable Battery Powered
Voltage dips IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/ 30 cycles Single phase: at 0°	i-Digits: Not applicable Charging base: 0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle and 70 % UT; 25/ 30 cycles Single phase: at 0°	Not applicable Battery Powered
Voltage interruptions IEC 61000-4-11	0 % UT; 250/300 cycle	i-Digits: Not applicable Charging Base: 0 % UT; 250/ 300 cycle	Not applicable Battery Powered
Power frequency (50/ 60 Hz) magnetic field IEC 61000-4-8	30 A/m 50/ 60Hz	i-Digits: Not applicable Charging Base: 30 A/m	Battery Powered

Guidance and manufacturer's declaration – electromagnetic immunity			
Össur Myoelectric Prosthetic Devices are intended for use in the electromagnetic environment specified below. The customer or the user of the Össur Myoelectric Prosthetic Devices should assure that each are used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted Disturbances Induced by RF Fields	IEC 61000-4-6	i-Digits; Not applicable Battery Powered No Cables >3m Charging Base; 3V 0.15 MHz – 80 MHz 6V in ISM and amateur bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	Portable and mobile RF communications equipment should be used no closer to any part of the Össur Myoelectric Prosthetic Devices including cables, than the recommended separation distance calculated from the equation appropriate to the frequency of the transmitter. Recommended separation distance d = 1.2 √P d = 1.2 √P 80 MHz to 800 MHz d = 2.3 √P 800 MHz to 2.7 GHz
Radiated RF IEC	IEC 61000-4-3	i-Digits; 12 V/m 26 MHz to 1GHz 10 V/m 1000 MHz to 2.7GHz 80% AM at 1kHz Charging Base; 10 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters as determined by an electromagnetic site survey ^a should be less than the compliance level in each frequency range ^b
Proximity magnetic fields	IEC 61000-4-39	Charging Base; 30kHz at 8 A/m As per table 1 of IEC 60601-1-2: 2014 +AMD 1: 2020	Interference may occur in the vicinity of equipment marked with the following symbol: 
IMMUNITY to proximity fields from RF wireless communications equipment	IEC 61000-4-3	Charging Base; As per table 9 of IEC 60601-1-2: 2014 +AMD 1: 2020	
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Össur Myoelectric Prosthetic Devices are used exceeds the applicable RF compliance level above, Össur Myoelectric Prosthetic Devices should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re- orienting or relocating the Össur Myoelectric Prosthetic Devices			

Test frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	27
450	430-470	GMRS 460, FRS 460	Pulse modulation 18 Hz	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	9
745				
780				
810	800-960	GSM 800/ 900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
870				
930				
1720	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, UMTS LTE Band 1, 3, 4, 25	Pulse modulation 217 Hz	28
1845				
1970				
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
5500				
5785				

If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3

- a) For some services, only the uplink frequencies are included.
- b) The carrier shall be modulated using a 50 % duty cycle square wave signal.
- c) As an alternative to FM modulation, the carrier may be pulse modulated using a 50 % duty cycle square wave signal at 18 Hz. While it does not represent actual modulation, it would be worst case

Recommended separation distance between portable and mobile RF communications equipment and the Össur Myoelectric Prosthetic Devices			
Össur Myoelectric Prosthetic Devices are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customers or the user of Össur Myoelectric Prosthetic Devices can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Össur Myoelectric Prosthetic Devices as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter in Watt	Separation distance according to frequency of transmitter in meters		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.7GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			

REPORT A SERIOUS INCIDENT

Important notice to users and/or patients established in Europe:

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

DESCRIPTION OF SYMBOLS

BF Applied Part

Manufacturer and year of manufacturing (YYYY)

Consult Instructions for use

Device contains electronic components and/or batteries that should not be disposed of in regular waste

IP22

Protected against solid foreign objects of 12.5mm diameter and greater and protected against splashing water.

Serial Number

For i-Digits™ devices:
The unique serial number for i-Digits™ devices is 2 letters followed by 6 digits.

European Conformity

Medical Device

Class II equipment – provides double isolation to protect against electric shock

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